 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
CONSTRUX® Mini PTC Spacer System

-  Click directory below for desired language
- | | | |
|----------------|--------------|------------|
| English | EN | 2-3 |
| Español | ES | 4-5 |
| 英文 | ZH-TW | 6-7 |

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

CONSTRUX® Mini PTC Spacer System

Description:

The CONSTRUX Mini PTC Spacer System is comprised of a variety of implants that have a PEEK core with integrated porous titanium end plates. CONSTRUX Mini PTC spacers are implanted in the cervical intervertebral disc space and are intended to facilitate vertebral fusion by stabilizing adjacent vertebrae, maintaining disc height and preventing the collapse of one vertebra onto another.

The CONSTRUX Mini PTC Spacer System is not intended to be used as a stand-alone device and must be used with supplemental fixation. The implants are used singly and are implanted using an anterior approach.

The CONSTRUX Mini PTC implants are provided sterile.

CONSTRUX Mini PTC implants are design to be used with CONSTRUX Mini PEEK Spacer System instrumentation. The CONSTRUX Mini PTC implants are not compatible with components or metal from any other manufacturer's system.

Indications for Use:

The CONSTRUX Mini PTC Spacer System is indicated for spinal fusion procedures at one or two contiguous levels within the cervical spine (C2-T1) in skeletally mature patients with degenerative disc disease (DDD). DDD is defined as neck pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies.

The CONSTRUX Mini PTC Spacer System is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and supplemental fixation system (i.e., an anterior cervical plate such as the Orthofix Hallmark Anterior Cervical Plate System).

Patients must have undergone a regimen of at least six weeks of non-operative treatment prior to being treated with the CONSTRUX Mini PTC Spacer System in the cervical spine.

Contraindications:

The CONSTRUX Mini PTC Spacer System, as with other orthopedic implants, is contraindicated for use in patients with:

1. Active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection.
2. Rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis which may prevent adequate fixation.
3. Conditions that may place excessive stresses on bone and implants such as severe obesity, pregnancy or degenerative diseases. The decision to use this system in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
4. Prior fusion at the level to be treated.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

Note: As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, permanent pain and/or deformity. In rare instances, some complications may be fatal.

Warnings and Precautions:

The surgeon should be aware of the following when using implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape, and design of the implant. No implant can be expected to withstand the unsupported stresses of full weight bearing. The size, shape and condition of human bones are also contributing factors to the success of the surgery.
2. **DO NOT USE DAMAGED IMPLANTS.** The correct handling of the implant is extremely important. Implants should not be bent, notched or scratched. These operations can produce defects in surface finish and may cause internal stress concentrations which may become the focal point for eventual failure of the device.
3. CONSTRUX Mini PTC implants are provided STERILE. Do not use if the package is opened or damaged or if the expiration date has passed.
4. **DO NOT re-sterilize** these implants as this could result in injury or require reoperation due to breakage.
5. **SINGLE USE ONLY** – CONSTRUX Mini PTC Spacer System implants are intended for Single Use Only. No surgical implants should be reused. Any implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns that may lead to fatigue failure. Reuse could result in injury or require reoperation due to breakage or infection.
6. **Non-Sterile;** the CONSTRUX Mini PEEK Spacer System instruments are provided non-sterile, and therefore, must be thoroughly cleaned and sterilized before initial use and after each subsequent use.
7. Postoperative care is important. The patient should be instructed in the limitations of the implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
8. Based on dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level and other patient conditions that may impact the performance of the intervertebral body fusion device.
9. The implantation of the intervertebral body fusion device should be performed only by experienced spinal surgeons with specific training in the use of the device because it is a technically demanding procedure presenting a risk of serious injury to the patient.

MRI Compatibility Information:

The CONSTRUX Mini PTC Spacer System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

Cleaning:

CONSTRUX Mini PTC Spacer System implants are provided STERILE. Do not use the implant if the package is opened or damaged or if the expiration date has passed. Please discard all open and unused implants. Do not re-sterilize an opened and unused implant.

The CONSTRUX Mini PEEK Spacer System instruments are compatible for use with the CONSTRUX Mini PTC Spacer System. All instruments must be thoroughly cleaned and sterilized after each use.

Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

None of the instruments require disassembly prior to cleaning.

From Point of Use:

Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Preparation for Cleaning:

1. Any instruments with moving parts (e.g., knobs, triggers and hinges) should be placed in the open position to allow access of the cleaning fluid to areas which are difficult to clean.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or a soft plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning. Use a soft plastic bristle brush or a pipe cleaner to remove soil from any inner lumens.

Manual Cleaning:

1. Upon completing the preparation for the cleaning procedure, prepare Vesphene® Ilse agent per the dilution recommended on the label directions – 1 ounce per gallon or 1 mL of Vesphene® Ilse to 128 mL of potable tap water per manufacturer’s Directions for Use label.
2. Bathe instruments in prepared room temperature solution.
3. Manually agitate instruments in Vesphene® Ilse solution for 15 minutes.
4. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of the interior lumen(s) to remove the soil.
5. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
6. Hang dry the device.
7. Visually inspect the instruments for visible soil.
8. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. Upon completing the preparation for the cleaning procedure, set up the washer/ disinfectant detergent dose at ½ ounce of Endozime AW Plus® per gallon of water or according to manufacturer’s recommendations.
2. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of any interior lumen(s) to remove the soil.
3. Place scrubbed instruments into the washer baskets.
4. Orient instruments in the automated washer’s carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
 - a. Pre-Wash 1: cold potable water, 2 minutes
 - b. Enzyme/Detergent treatment:
 1. Spray with enzyme/detergent, 20 seconds
 2. Soak, 1 minute
 3. Rinse cold potable water, 15 seconds
 4. Rinse cold potable water, 15 seconds
 - c. Wash ≥ 65°C, 2 minutes using a Endozime AW Plus®
 - d. Rinse 1: hot potable water, 15 seconds
 - e. Rinse 2: hot potable water, 15 seconds
 - f. Rinse 3: hot potable water, 15 seconds
 - g. Rinse 4: hot potable water, 15 seconds
 - h. Thermal rinse ≥ 93°C, 1 minute
 - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
 - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
 - k. Dry at 115°C, 7 minutes
6. Inspect the instruments for visible soil.
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix any instruments that are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

Sterilization:

CONSTRUX Mini PTC Spacer System implants are provided STERILE. They are sterilized using gamma irradiation sterilization. Do not re-sterilize.

Sterilization in Orthofix Cases with Blue Wrap:

The CONSTRUX Mini PEEK Spacer System instruments are supplied NON-STERILE and are compatible for use with the CONSTRUX Mini PTC Spacer System. Prior to use, all instruments should be placed in the appropriate Orthofix case which will be wrapped in a FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Preconditioning: 4 pulses
Drying time: 30 minutes		Exposure time: 4 minutes
Double wrapped		Drying time: 30 minutes
		Double wrapped

Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer’s instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix case (recommended 23¼” long x 11¼” wide container). The following sterilization cycle has been validated:

Method: Steam
Cycle: Prevac
Temperature: 270°F (132°C)
Preconditioning: 4 pulses
Exposure time: 4 minutes
Drying time: 30 minutes

Note: Rigid sterilization containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed per ANSI/AAMI ST79 Comprehensive guide to steam sterilization and sterility assurance in health care facilities. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The CONSTRUX Mini Spacer System instruments are provided in modular cases specifically intended to contain and organize the system’s components. The system’s instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness, and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.



Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the CONSTRUX Mini PTC Spacer System are not formulated with and do not contain natural rubber. The term “natural rubber” includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Rx Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use	 Use By Date
	Orthofix.com/IFU	 Manufacturer
	Single Use Only Do Not Reuse	 Authorized Representative
	Catalogue Number	 Do Not Resterilize
	Provided Non-Sterile	 Serial Number
	Sterilized Using Irradiation	 Lot Number

Orthofix Inc.
 3451 Plano Parkway
 Lewisville, Texas 75056-9453 U.S.A.
 1-214-937-3199
 1-888-298-5700
 www.orthofix.com

Australian Sponsor
 Emergo Australia
 Level 20, Tower II
 Darling Park
 201 Sussex Street
 Sydney, NSW 2000
 Australia

EC/REP Medical Device Safety Service
 (MDSS)
 Schiffgraben 41
 30175 Hannover
 Germany
 +49 511 6262 8630
 www.mdss.com

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