 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
InSWing™ Interspinous Spacer

 Click directory below for desired language

English	EN	2-3
Español	ES	4-5
Français	FR	6-7
Italiano	IT	8-9

Orthofix Inc.
 3451 Plano Parkway
 Lewisville, Texas 75056-9453 U.S.A.
 1-214-937-3199
 1-888-298-5700
 www.orthofix.com

Australian Sponsor
 Emergo Australia
 Level 20, Tower II
 Darling Park
 201 Sussex Street
 Sydney, NSW 2000
 Australia

EC REP Medical Device Safety Service
 (MDSS)
 Schiffgraben 41
 30175 Hannover
 Germany
 +49 511 6262 8630
 www.mdss.com

English EN

Device System Name:

InSWing™ Interspinous Spacer

Description:

The InSWing Interspinous Spacer is an implantable spinal construct that is manufactured from PEEK and is designed to fit between the spinous processes of the vertebrae. InSWing is implanted through a unilateral approach. The self-deploying and self-locking wings are designed to prevent the implant from pulling out. The implant limits extension and the flexion-limiting band (optional) is intended to improve the stability of the treated intervertebral lumbar segment.

The system is composed of the following components, each provided in an individual sterile package:

- InSWing Spacer implant 8 mm, Catalogue No. 17-0008
- InSWing Spacer implant 10 mm, Catalogue No. 17-0010
- InSWing Spacer implant 12 mm, Catalogue No. 17-0012
- InSWing Spacer implant 14 mm, Catalogue No. 17-0014
- InSWing Spacer implant 16 mm, Catalogue No. 17-0016

The following components are provided sterile and are packaged together as Catalogue No. 17-0001: InSWing band (qty. of 1), and clip (qty. of 2).

The InSWing spacers are made out of PEEK (polyetheretherketone) as described by ASTM F2026, with titanium alloy (Ti-6Al-4V) markers as described by ASTM F136. The InSWing band is made out of 100% PET (polyethylene terephthalate). The InSWing clip is made out of titanium alloy (Ti-6Al-4V) as described by ASTM F136.

The InSWing Interspinous Spacer is intended to be placed between the spinous processes of the vertebrae via a unilateral approach.

Indications for Use:

The InSWing Interspinous Spacer is indicated to relieve leg and back pain due to symptomatic Lumbar Spinal Stenosis (LSS) from levels L1/L2 to L4/L5 with or without neurogenic intermittent claudication. The device can be used in combination with decompressive procedures providing that the spinous processes remained intact.

Contraindications:

Contraindications include, but are not limited to:

1. Morbid obesity.
2. Mental illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia/osteoporosis.
7. Patients unwilling or unable to follow post-operative care instructions.
8. Active local and systemic inflammation conditions.
9. Patient allergy to implant material or to its elements.
10. Patient's poor general condition that has a high risk of endangering health condition and/or loss of life as a result of surgical procedure.
11. Paget's disease.
12. Kyphosis requiring surgical correction.
13. Any circumstances not listed under the heading Indications for Use.
14. Other general or patient related contraindications against surgical implant treatment which may preclude potential benefit from implantation.

Potential Adverse Events:

All of the possible adverse events associated with spinal surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Neurological injury.
2. Vascular or visceral injury.
3. Fracture of the spinous process.
4. Early or late loosening of any or all of the components.
5. Incorrect positioning of the implant.
6. Disassembly, breakage and/or bending of any or all components.
7. Foreign body (allergic) reaction to implants, debris, corrosion products, including metallosis, straining, tumor formation, and/or auto-immune disease.
8. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
9. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.

10. Infection.
11. Pain, discomfort, or abnormal sensations due to the presence of the device.
12. Chronic inflammation.
13. Hemorrhage.
14. Cessation of any potential growth of the operated portion of the spine.
15. Death.
16. Potential risks identified with the use of the device system may require additional surgery.

Warnings and Precautions:

1. The safety and effectiveness of the interspinous spacer have been established only for spinal conditions with symptomatic LSS of the L1/L2-L4/L5 levels.
2. Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, neurological injury, and vascular or visceral injury.
3. Single use only.
4. Non-sterile: instruments are sold non-sterile, and therefore must be sterilized before use.
5. DO NOT REUSE IMPLANTS. Discard used, damaged, or otherwise suspect implants.
6. Do not use the InSWing implant past the expiration date noted on the package.
7. Do not use damaged products or opened packages.
8. Do not re-sterilize implants.
9. The implant should be placed as close as possible to the lamina without removing the supraspinous ligament. The interspinous ligament should be properly removed, as per the InSWing operative technique.
10. The implantation of the InSWing Interspinous Spacer System should be performed only by experienced spinal surgeons with specific training in the use of this implant because it is a technically demanding procedure presenting a risk of serious injury to the patient.
11. Do not use the components of this system with implants of other material composition or components from different manufacturers unless specifically stated.
12. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not attempt to re-sterilize single-use implants.

MRI Compatibility Information:

The InSWing Interspinous Spacer has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

Cleaning:

All non-sterile instruments must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

The InSWing Interspinous Spacer implants are provided sterile and do not require cleaning.

Sterilization:

The InSWing Interspinous Spacer instruments (medical-grade stainless steel and silicone handles) should be steam sterilized by the hospital using the following recommended cycles when utilizing a sterilization wrap:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

The InSWing Interspinous Spacer implants are sterilized using gamma irradiation sterilization. Do not re-sterilize.

Physician Information:

Preoperative:

1. Only patients that meet the criteria described in the indications should be selected.
2. Patient conditions and/or predispositions such as those addressed in the aforementioned contraindications should be avoided.
3. Correct selection of the implant is extremely important.
4. Use care in handling and storage of implant components. Inspection should be made to determine if components have been damaged during storage or previous procedures.
5. An adequate inventory of implant sizes should be available at the time of surgery.
6. Certain special surgical instruments are required to perform this surgery. Review of the use and handling of these instruments is very important.

Intraoperative:

1. A unilateral surgical approach is the recommended access to the interspinous space in order to minimize trauma to the soft tissue.
2. Use great care to insure that the spinous process scraper is used to clean the remains of the interspinous ligament from spinous process.
3. Do not damage the supraspinous ligament. Be sure to detach muscle on the contralateral side with the scraper to allow the wings to properly deploy
4. Correct sizing is imperative for good clinical outcome. When using the sizing trial or distractor, be careful not to over distract. Too much distraction can induce kyphosis.
5. During implantation of the InSWing, do not use excessive force as it may cause a fracture of the spinous process.
6. Before closing, use fluoroscopy to verify the final implant position. Complete the surgery using standard posterior lumbar closure procedures.

Postoperative:

1. Detailed instructions on the use and limitations of the implant should be given to the patient. The patient must be made aware of the limitations of the implant.
2. Periodic X-rays for at least the first year postoperatively are recommended for close comparison with postoperative conditions to detect any evidence of changes in position, bending or cracking of components.
3. Surgical implants must never be reused. Any retrieved devices should never be reused in another surgical procedure. The retrieved parts should be handled and disposed of in such a manner as to ensure that reuse is not possible.

Removal of Implant:

In certain circumstances, removal may be necessary. These reasons include, but are not limited to:

1. Dismantling or breakage of implant.
2. Fracture of spinous process with pain.
3. Infection.
4. Surgery involving spinous processes.

The recommended surgical steps for removal are as follows:

1. Cut the band and knot (if present).
2. Remove the band and clip used to secure the original knot.
3. Make a small opening in the intact fascia on the side of the opening wings.
4. Insert a raspatory (periosteal elevator) along the spinous process and close (lift) one of the wings of the implant while pulling the implant out. This can also be done with a finger.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The InSWing Interspinous Spacer implants and instruments should be stored in a cool and dry place.

The InSWing Interspinous Spacer instruments are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

The InSWing Interspinous Spacer implants are supplied sterile and packaged in a dual barrier system to facilitate introduction into the sterile field in the surgical theater. The devices are considered sterile, unless the sterile expiration date is elapsed or the package is opened or damaged, e.g., loss of seal integrity or is not in the original packaging.

Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.



Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the InSWing Interspinous Spacer System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: The InSWing Interspinous Spacer has not been cleared by the US Food and Drug Administration and is not available in the United States.

	See Instructions for Use		Use By Date
	Orthofix.com/IFU		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Catalogue Number		Do Not Resterilize
	Lot Number		Serial Number
	Sterilized Using Irradiation		

Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

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