<p>ORTHOFIX® FIXATION SYSTEM MEDICAL DEVICES</p> <p>3</p>	EN	<p>MEDICINSKA ANORDNINGAR I ORTHOFIX®FIXERINGSSYSTEM</p> <p>66</p>	SV
<p>DISPOSITIVI MEDICALI DEL SISTEMA DI FISSAZIONE ORTHOFIX®</p> <p>9</p>	IT	<p>ΟΔΗΓΙΕΣ ΓΙΑ ΤΗΝ ΑΣΦΑΛΗ ΕΠΙΞΕΡΓΑΣΙΑ ΤΩΝ ΙΑΤΡΙΚΩΝ ΣΥΣΚΕΥΩΝ ΤΟΥ ΣΥΣΤΗΜΑΤΟΣ ΣΤΑΘΕΡΟΠΟΙΗΣΗΣ ORTHOFIX</p> <p>72</p>	EL
<p>DISPOSITIFS MÉDICAUX DU SYSTÈME DE FIXATION ORTHOFIX®</p> <p>16</p>	FR	<p>医療用固定システムを使用する為 の安全な手順について</p> <p>79</p>	JP
<p>MEDIZINGERÄTE DES ORTHOFIX® FIXATIONSSYSTEMS</p> <p>23</p>	DE	<p>ORTHOFIX® 固定系统医疗器械</p> <p>85</p>	CN
<p>DISPOSITIVOS MEDICOS DEL SISTEMA DE FIJACION ORTHOFIX®</p> <p>30</p>	ES	<p>MEDICÍNSKÁ ZAŘIZENÍ A NÁSTROJE FIXAČNÍHO SYSTÉMU ORTHOFIX®</p> <p>91</p>	CS
<p>ORTHOFIX® MEDICINSK UDSTYR TIL FIKSATION</p> <p>36</p>	DA	<p>ORTHOFIX® PRZYRZĄDY MEDYCZNE – SYSTEM STABILIZACJI ZEWNĘTRZNEJ</p> <p>98</p>	PL
<p>OHJEET ORTHOFIX® - KIINNITYSJÄRJESTELMÄÄN KUULUVIEN LÄÄKETIETEELLISTEN LAITTEIDEN TURVALLISTA KÄSITTELYÄ VARTEN</p> <p>42</p>	FI	<p>MEDICINSKIH PRIPOMOČKOV ORTHOFIX® FIKSACIJSKEGA SISTEMA</p> <p>105</p>	SL
<p>ORTHOFIX® MEDISINSK UTSTYR FOR FIKSERINGSSYSTEMER</p> <p>48</p>	NO	<p>ORTHOFIX® 고정 기구 의료 기기</p> <p>111</p>	KR
<p>STUKKEN VAN HET ORTHOFIX® FIXATIESYSTEEM</p> <p>54</p>	NL		
<p>DISPOSITIVOS MÉDICOS DO SISTEMA DE FIXAÇÃO ORTHOFIX®</p> <p>60</p>	PT	<p>ORTHOFIX Srl Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy Telephone 0039-0456719000 - Fax 0039-0456719380</p>  0123	PQ ISP M 03/12

# Instructions for the Safe Processing of the ORTHOFIX® FIXATION SYSTEM MEDICAL DEVICES

<b>Manufacturer name</b>	<b>ORTHOFIX SRL</b> via Delle Nazioni 9 37012 Bussolengo (VR) Italy Tel 0039 (0)45 6719000 fax 0039(0)45 6719380
<b>Method code</b>	<b>OF1: processing of the Orthofix® Fixation System medical devices</b>
<b>Symbol</b>	 Attention, see instructions for use PQ ISP
<b>Device(s)</b>	All the Orthofix, Fixation System products NOT SUPPLIED "STERILE"
<b>WARNINGS</b>	Aluminium based instruments are damaged by alkaline (pH>7) detergents and solutions. Anodised coating is damaged by detergents with free halogen ions or sodium hydroxide. Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
<b>Limitations and restrictions on reprocessing</b>	PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED. Repeated reprocessing has minimal effect on reusable instruments. End of life is normally determined by wear and damage due to use.

## INSTRUCTIONS FOR PROCESSING NEW DEVICES SUPPLIED "NON-STERILE" PRIOR TO THEIR FIRST USE

<b>General</b>	<ul style="list-style-type: none"> <li>• Unless supplied sterile, all Orthofix medical devices must be sterilized prior to surgical use.</li> <li>• A new product means any device taken out of its original Orthofix packaging.</li> </ul>
<b>Recommended decontamination process</b>	The following sequence of processes is recommended for rendering the Orthofix Fixation System devices which are supplied non-sterile safe for their <b>first clinical use</b> : 1. CLEANING 2. STERILIZATION

### Preparation for cleaning of new products

- When applicable, the device should be disassembled completely.
- Wherever possible, all parts of disassembled devices should be kept together in one container.
- For disassembly/reassembly the specific instrumentation should be used as described in detail in the Operative Manuals.

### STEP 1: Cleaning of new products

- REMOVE products from their original packaging. All equipment should be carefully examined prior to use to assure proper working condition.
- CLEAN with a woven-non woven tissue soaked using a solution of 70% medical grade alcohol and 30% distilled water or with compatible detergent.  
**Detergents with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.**
- RINSE with sterile distilled water.

### STEP 2: Drying

- HAND-DRY carefully, using absorbent, non-shedding cloth or industrial hot hair dryer, or place into drying cabinet.

### STEP 3: Sterilization of new products

- Prior to surgical use, new products should be cleaned as described at Step 1.
- PACKAGING: Where products are to be packaged to maintain their sterility after sterilization and to prevent damage of the instrument prior to use, an appropriate medical grade packaging material should be used. Ensure that the pack is large enough to contain the instruments without stressing the seals.
- STERILIZATION CONTAINERS: Instruments may be loaded into a dedicated instrument tray, or general-purpose sterilization tray. Ensure that cutting edges are protected and do not exceed the recommended content or maximum weight indicated by manufacturer.
- PRECAUTIONS: Fixators can be sterilized in the assembled state as long as ball-joints, central body locking nut and clamp locking screws are left untightened. If any joints are tightened, they may sustain damage from thermal expansion during the sterilization process.
- STERILIZATION: Sterilize by steam autoclaving, utilizing a pre-vacuum cycle.  
*Orthofix recommends the following cycle:  
Steam autoclave 132-135°C (270-275°F), minimum holding time 10 minutes.*
- Any other **validated** pre-vacuum autoclave cycle may be used in alternative.

## INSTRUCTIONS FOR REPROCESSING OF REUSABLE DEVICES FOR SUBSEQUENT RE-USE

### General

- Unless supplied sterile, all Orthofix devices must be sterilized prior to surgical use.
- PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED.
- Repeated reprocessing has minimal effect on reusable instruments. End of life is normally determined by wear and damage due to use.

### Recommended decontamination process

The following sequence of decontamination processes is recommended for **reprocessing re-usable Orthofix Fixation System medical devices** and rendering them safe for subsequent clinical use.

1. CLEANING
2. DISINFECTION
3. STERILIZATION

### Preparation at the point of use of used devices

- All **used** surgical instruments should be regarded as **contaminated**. Their handling, collection and transportation should be strictly controlled to minimise any possible risks to patients, personnel and any area of the healthcare facility.
- It is recommended that instruments are reprocessed as soon as is reasonably practicable following use.
- Ensure that items are securely and safely packaged during transport to the decontamination area.

## Preparation for cleaning of used devices

- Used fixators should be disassembled completely, and all cams and bushes discarded since they are not re-usable.
- All components should be inspected, since damage to the surface of metal components can reduce the strength and fatigue resistance and may lead to corrosion.
- Wherever possible, all parts of disassembled devices should be kept together in one container.
- For disassembly/reassembly, the specific instrumentation should be used as described in detail in the Operative Manuals.

## STEP 1: Cleaning of used products

- Cleaning is an essential pre-requisite to ensure effective disinfection or sterilization.
- The preferred method of decontamination of used devices is mechanical cleaning followed by disinfection.
- Where an automated washer-disinfector is not available, manual cleaning may be used, followed by disinfection.

### MANUAL CLEANING (BY IMMERSION)

#### *Equipment required:*

- a sink (not hand wash basin) or receptacle which will hold enough detergent so that the item of equipment to be cleaned can be fully immersed;
- a detergent solution. **Orthofix recommends use of a 0,3% enzymatic detergent solution, immersion for 30 minutes at 40°C (104°F). Detergents with fluoride, chloride, bromide, iodide, or hydroxyl ions (free halogen ions or sodium hydroxide) MUST NOT be used;**
- a receptacle to contain rinse water;
- a drainage surface;
- a clean, disposable, absorbent, non-shedding cloth or mechanical drying facility (drying cabinet or industrial hot air dryer);
- a brush and jet washer.

#### *Procedure:*

- I. Ensure that the cleaning receptacle is clean and dry.
- II. Wearing protective equipment, fill the receptacle with sufficient water/detergent solution.
- III. Carefully immerse all components in the solution in order to displace trapped air; it is important to ensure that the cleaning solution reached all surfaces, including those of devices that have holes or recesses or are cannulated.
- IV. Brush, wipe, agitate, irrigate, jetwash or hand spray the item to dislodge and remove all visible dirt, ensuring that the action is performed beneath the surface of the solution.
- V. Remove the items from the solution and drain.
- VI. Remove any residue with a brush in running water.
- VII. Soak in sterile distilled water to remove traces of hard water.
- VIII. Remove item from rinse water and drain.
- IX. Carefully hand-dry using absorbent, non-shedding cloth or an industrial hot air dryer, or place in a drying cabinet.
- X. Complete the necessary documentation.
- XI. Proceed with disinfection.

### MECHANICAL CLEANING (USING AN AUTOMATED WASHER)

- Mechanical cleaning followed by disinfection is the preferred method of decontamination of used devices.
- If a washer-disinfector is used, it must have a validated cycle.
- Ensure that the washer-disinfector and all services are operational.
- Select and start a cycle according to the recommendations of the washer manufacturer
- **Detergents with fluoride, chloride, bromide, iodide, or hydroxyl ions (free halogen ions or sodium hydroxide) MUST NOT be used.**
- Proceed with disinfection.

**STEP 2:**  
**Disinfection**

**DISINFECTION PROCEDURE (MANUAL)**

*Equipment required:*

- a sink (not hand wash basin) or receptacle which will hold enough disinfectant so that the item of equipment to be cleaned can be fully immersed;
  - a compatible water/disinfectant solution at dilution and temperature recommended by its producer. **Disinfectants with fluoride, chloride, bromide, iodide, or hydroxyl ions [free halogen ions or sodium hydroxide] MUST NOT be used; Orthofix recommends use of 3% hydrogen peroxide, with immersion, for 3 hours at room temperature;**
  - a receptacle to contain rinse water;
  - a drainage surface;
  - a clean, disposable, absorbent, non-shedding cloth or mechanical drying facility (drying cabinet or industrial hot air dryer);
  - a brush and jet washer.
- 

*Procedure:*

- I. Wearing protective equipment, fill the receptacle with sufficient disinfectant solution to ensure complete immersion of the item.
  - II. Carefully immerse all components in the solution in order to displace trapped air; it is important to ensure that the cleaning solution reached all surfaces, including those of cannulated devices.
  - III. Leave the items for the time required (3 hours with 3% Hydrogen Peroxide).
  - IV. Remove the items from the solution and drain.
  - V. Remove any residue with a brush in running water.
  - VI. Soak in sterile distilled water to remove traces of hard water.
  - VII. Remove item from rinse water and drain.
  - VIII. Carefully hand-dry using absorbent, non-shedding cloth or industrial hot hair dryer, or place in a drying cabinet.
  - IX. Complete the necessary documentation.
  - X. Proceed with sterilization.
- 

**DISINFECTION PROCEDURE (AUTOMATIC)**

*Equipment required:*

- I. A thermal washer-disinfector, cabinet or continuous process type.
  - II. A sufficient number of racks for stacking items to be processed.
  - III. A **compatible** disinfectant and rinse aid. **Disinfectants with fluoride, chloride, bromide, iodide, or hydroxyl ions [free halogen ions or sodium hydroxide] must not be used; Orthofix recommends use of 0,5% phenolic disinfectant solution at 80°C (176°F). The validated cycle time is 80 minutes.**
  - IV. A drainage surface.
  - V. A clean, disposable, absorbent, non-shedding cloth or mechanical drying facility.
- 

*Procedure:*

- I. Ensure the washer-disinfector and all services are operational.
- II. Wear protective equipment, load the rack/machine ensuring that the loading configuration does not impede the cleansing process.
- III. Select and start a cycle according to the recommendation of the washer manufacturer. On completion on the cycle, ensure that all stages and parameters have been achieved; remove the load and visually check and inspect the cleanliness of the item, drain off excessive water and dry if necessary.
- IV. Complete the documentation.
- V. Proceed with sterilization.

## DISINFECTION PROCEDURE (ULTRASOUND)

### Equipment required:

- I. An ultrasonic washer with lid which will hold enough liquid so that the items of equipment to be cleaned can be fully immersed.
- II. A sufficient number of supporting racks or trays for stacking items to be processed.
- III. A timing device.
- IV. A compatible water-detergent solution at dilution and temperature, recommended by manufacturer.
- V. A clean, disposable, absorbent, non-shedding cloth or mechanical drying facility.

---

### Procedure:

- I. Ensure the ultrasonic washer is clean and dry prior to use.
- II. Wear protective equipment, fill the fluid reservoir with sufficient water/disinfectant to ensure complete immersion of items. **Disinfectants with fluoride, chloride, bromide, iodide, or hydroxyl ions [free halogen ions or sodium hydroxide] must not be used; Orthofix recommends use of 0.5% phenolic disinfectant solution immersion at 50°C (122°F) for 15 minutes (ultrasound frequency 50/60 Hz).**
- III. Switch on and leave for required time to degas the water.
- IV. Remove lid and carefully immerse the item in the fluid ensuring that any air contained within the item is displaced. Irrigate cannulated devices.
- V. Re-place the lid and leave for the time recommended (15 minutes).
- VI. Switch off, lift the lid, remove the item and drain before transferring to a clean-rinse receptacle.
- VII. Rinse thoroughly with clean water, ensuring irrigation of lumen devices, and drain.
- VIII. Carefully hand-dry using absorbent, non-shedding cloth, industrial hot air dryer or place in a drying cabinet.
- IX. Complete the documentation.
- X. Proceed with sterilization.

---

## STEP 3: Drying

Carefully hand-dry using absorbent, non-shedding cloth or industrial hot air dryer, or place in a drying cabinet.

## STEP 4: Inspection, maintenance and testing

- All instruments and product components should be visually inspected for cleanness and any signs of deterioration that may cause failure in use (such as cracks or damage to surfaces) and functions tested before being sterilized (see detailed Operative Technique Manuals and Instructions for use). Particular attention should be given to:
  - Cannulated devices (NB: cannulated drill-bits are single-patient use only)
  - Cutting edges: Discard blunt or damaged instruments
  - Hinged instruments: check for smooth movement of hinges without excessive "play".
  - Locking mechanisms should be checked for action.
- If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate all parts, except for cam, bush and ball-joint coupling with lubrication oil for medical applications whenever required (see detailed Operative Technique Manuals).
- Final locking of the ball-joints of the 10000, 30000/31000 or 90000 range of fixators is performed with a torque wrench, which must be turned in a clockwise direction only. A click indicates that the correct torque has been applied. Any attempt to unlock the cam or any screw using the torque wrench will damage its gearing. The torque wrench is pre-set at a specific value, which is 15 Nm±0.5 for the 30000/31000 Range (30025) and 27 Nm±1 for the 10000 and 90000 Ranges (10025). This value should be checked at least every two years or any time the instrument becomes damaged, by returning it to the local authorized representative.  
The cams and bushes in the ball joints of the 10000, 30000, 31000 and 90000 range of fixators must be replaced after every use.

The use of a torque wrench is not required for final locking of the ball joints in the XCaliber External Fixator. Final locking is achieved with the Allen Wrench. The cams can be locked from either side of the clamp. They should be turned towards the thicker section of the coloured insert until tightly closed, and the cam is at least 50% of the way across the recess.

#### STEP 5: Sterilization

- **PACKAGING:** Where products are to be packaged to maintain their sterility after sterilization and to prevent damage of the instrument prior to use, an appropriate medical grade packaging material should be used. The pack should be large enough to contain the instruments without stressing the seals.
- **STERILIZATION CONTAINERS:** Instruments may be loaded into a dedicated (Orthofix) instrument tray, or general-purpose sterilization tray. Cutting edges should be protected and the recommended content or maximum weight not exceeded as indicated by manufacturer.
- **PRECAUTIONS:** Fixators can be sterilized in the assembled state as long as ball-joints, central body locking nut and clamp locking screws are left *untightened*. If any joints are tightened, they may sustain damage from thermal expansion during the sterilization process.
- **STERILIZATION:** Sterilize by steam autoclaving, utilizing a pre-vacuum cycle.  
*Orthofix recommends the following cycle:*  
*Steam autoclave 132-135°C (270-275°F), minimum holding time 10 minutes.*
- Any other **validated** pre-vacuum autoclave cycle may be used in alternative.
- Processed items should be stored in clean and secure stores to avoid damage or tampering.

**Disclaimer: "The instructions provided above have been validated by Orthofix as being a true description of the preparation of a device for first clinical use or for re-use of multiple use devices. It remains the responsibility of the reprocessor to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes should be adequately recorded. Likewise any deviation by the reprocessor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences and should also be appropriately recorded".**

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All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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