 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**FORZA® PTC Spacer System**

 Click directory below for desired language

<b>English</b> 	<b>2-4</b>
<b>Español</b> 	<b>5-7</b>
<b>Deutsch</b> 	<b>8-10</b>

 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

## English **EN**

### Device System Name:

## **FORZA<sup>®</sup> PTC Spacer System**

### Description:

The FORZA PTC Spacer System is comprised of a variety of implants that have a PEEK (OPTIMA LT1) core with integrated porous titanium (Ti-6Al-4V) end plates as well as a tantalum marker that acts as a visual aid for the surgeon in determining the location of the implant both intraoperatively and postoperatively.

FORZA PTC Spacer System implants are offered in two geometric shapes – straight and curved, and offered in parallel and lordotic profiles to restore the natural curvature of the spine. The implants can be used in single placement or in pairs. Both the curved and straight implants feature a bulleted nose for ease of insertion and a roughened surface on both the inferior and superior faces of the implant to provide increased stability and help prevent anterior/posterior movement of the device.

The FORZA PTC Spacer System is intended for intervertebral body fusion to aid in the surgical correction of the spine and are implanted using a posterior approach.

The FORZA PTC Spacer System is not intended to be used as a standalone device. The FORZA PTC Spacer System must be used with a supplemental fixation system.

The FORZA PTC Spacer System implants are provided sterile.

FORZA PTC implants are designed for use with FORZA PEEK Spacer System instrumentation. The FORZA PTC spacers are not compatible with components or metal from any other manufacturer's system.

### Indications for Use:

The FORZA PTC Spacer System is indicated for spinal fusion procedures in skeletally mature patients with degenerative disc disease (DDD) at one or two contiguous levels in the lumbar spine (L2-S1). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies. DDD patients may also have up to Grade 1 spondylolisthesis at the involved levels. These patients may have had a previous non-fusion surgery at the involved levels.

The FORZA PTC Spacer System is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and supplemental fixation, e.g., the Firebird Spinal Fixation System.

Patients must have undergone a regimen of at least six months of non-operative treatment prior to being treated with the FORZA PTC Spacer System.

### Contraindications:

The FORZA PTC Spacer System, as with other orthopedic implants, is contraindicated for use in patients with:

1. Active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection.
2. Rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis which may prevent adequate fixation.
3. Conditions that may place excessive stresses on bone and implants, such as severe obesity, pregnancy or degenerative diseases. The decision to use this system in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
4. Prior fusion at the level to be treated.
5. Any circumstances not listed under the heading Indications for Use.

### Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

**Note:** As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, permanent pain and/or deformity. In rare instances, some complications may be fatal.

### Warnings and Precautions:

The surgeon should be aware of the following when using implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. No implant can be expected to withstand the unsupported stresses of full weight bearing. The size, shape and condition of human bones are also contributing factors to the success of the surgery.
2. **DO NOT USE DAMAGED IMPLANTS.** The correct handling of the implant is extremely important. Implants should not be bent, notched or scratched. These operations can produce defects in surface finish and may cause internal stress concentrations which may become the focal point for eventual failure of the device.
3. FORZA PTC implants are provided **STERILE**. Do not use if the package is opened or damaged or if the expiration date has passed.
4. **DO NOT** re-sterilize these implants as this could result in injury or require reoperation due to breakage.
5. **SINGLE USE ONLY** – FORZA PTC Spacer System implants are intended for **SINGLE USE ONLY**. No surgical implants should be reused. Any implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns that may lead to fatigue failure. Reuse could result in injury or require reoperation due to breakage or infection.
6. **Non-Sterile;** the FORZA PEEK Spacer System instruments are provided non-sterile, and therefore, must be thoroughly cleaned and sterilized before initial use and after each subsequent use.
7. Postoperative care is important. The patient should be instructed in the limitations of the implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
8. Based on dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level and other patient conditions that may impact the performance of the intervertebral body fusion device.
9. The implantation of the intervertebral body fusion device should be performed only by experienced spinal surgeons with specific training in the use of the device because it is a technically demanding procedure presenting a risk of serious injury to the patient.
10. Patients with previous surgery at the levels to be treated may have different clinical outcomes compared with those without a previous surgery.

### MRI Compatibility Information:

The FORZA PTC Spacer System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

### Cleaning:

FORZA PTC Spacer System implants are provided **STERILE**. Do not use the implant if the package is opened or damaged or if the expiration date has passed. Please discard all opened and unused implants. Do not re-sterilize an opened and unused implant.

The FORZA PEEK Spacer System instruments are compatible for use with the FORZA PTC Spacer System. All instruments must be thoroughly cleaned and sterilized after each use.

Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

FORZA modular Implant Inserters require disassembly prior to cleaning. None of the other FORZA Spacer System instruments require disassembly prior to cleaning.

### From Point of Use:

Whenever possible, do not allow the blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

### Preparation for Cleaning:

1. All instruments with moving parts (e.g., knobs, triggers and hinges) should be placed in the open position to allow access of the cleaning fluid to areas which are difficult to clean.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or a soft plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning. Use a soft plastic bristle brush or a pipe cleaner to remove soil from any inner lumens.

### Instructions for Disassembly:

1. Modular Implant Inserters:  
There are three parts to each modular Implant Inserter: 1. Inserter Shaft, 2. Instrument Body and 3. Threaded Insert. The Universal Implant Inserter Handle may also be attached. The modular Implant Inserters require disassembly prior to cleaning using the following steps:
  - a. Grip the instrument at the gold-colored knob and Inserter Shaft with one hand.
  - b. With the other hand, hold the black silicone handle and pull the plunger away from the Instrument Body.
  - c. While still pulling the plunger, rotate the handle clockwise until disengaged from the Instrument Body. Set aside the Universal Implant Inserter Handle.
  - d. Hold the assembly level and grip the Instrument Body at the gold-colored knob with the pronged tip pointed away from you.
  - e. Turn the ribbed collar on the Inserter Shaft and firmly rotate it clockwise to disengage it from the Instrument Body. Completely unthread the ribbed collar until it slides down and stops on colored ring.
  - f. Keep the ribbed collar away from the threaded side near colored ring and firmly tug the Inserter Shaft with the other hand. Use caution to not drop the Threaded Insert out of the Inserter Shaft.
  - g. Slowly tilt the Inserter Shaft so the prongs are higher than the ribbed collar to remove the Threaded Insert. If the Threaded Insert does not disengage, gently tap the larger opening onto the palm of your hand to remove it.

### Manual Cleaning:

1. Upon completing the preparation for the cleaning procedure, prepare the Vesphene® Ise agent per the dilution recommended on the label directions – 1 ounce per gallon or 1 mL of Vesphene® Ise to each 128 mL of potable tap water per the manufacturer's Directions for Use label.
2. Bathe instruments in prepared room temperature solution.
3. Manually agitate instruments in Vesphene® Ise solution for 15 minutes.
4. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of any interior lumen(s) to remove the soil.
5. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
6. Hang dry the device.
7. Inspect the instruments for visible soil.
8. If visible soil is noted, repeat the steps listed above.

### Automated Cleaning:

1. Upon completing the preparation for the cleaning procedure, set up the washer/disinfectant detergent dose at ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer's recommendations.
2. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of any interior lumen(s) to remove the soil.
3. Place scrubbed instruments into the washer baskets.
4. Orient instruments in the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
  - a. Pre-Wash 1: cold potable water, 2 minutes
  - b. Enzyme/Detergent treatment:
    1. Spray with enzyme/detergent, 20 seconds
    2. Soak, 1 minute
    3. Rinse cold potable water, 15 seconds
    4. Rinse cold potable water, 15 seconds
  - c. Wash ≥ 65°C, 2 minutes using Endozime AW Plus®
  - d. Rinse 1: hot potable water, 15 seconds
  - e. Rinse 2: hot potable water, 15 seconds
  - f. Rinse 3: hot potable water, 15 seconds
  - g. Rinse 4: hot potable water, 15 seconds
  - h. Thermal rinse ≥ 93°C, 1 minute
  - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
  - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
  - k. Dry at 115°C, 7 minutes
6. Inspect the instruments for visible soil.
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

**Note:** Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

**Note:** Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix any instruments that are broken, discolored, corroded, have cracked components, pits, gouges, or are otherwise found defective. Do not use defective instruments.

### Instructions for Assembly and Inspection:

1. Modular Implant Inserters:  
There are three parts to each modular Implant Inserter: 1. Inserter Shaft, 2. Instrument Body and 3. Threaded Insert. The modular Implant Inserters require assembly after cleaning using the following steps:
  - a. Hold the Inserter Shaft so that the prongs are pointed toward the floor. Carefully insert the Threaded Insert, threaded tip first, into the large opening on the Inserter Shaft.
  - b. Keep the prongs pointed toward the floor with the ribbed collar resting on the colored ring in one hand, and insert the Instrument Body with the gold-colored knob into the large opening on Inserter Shaft.
  - c. Turn the Instrument Body until it fully seats in the in the square connection. There should be no gap between the Instrument Body and Inserter Shaft.
  - d. With the prongs still pointed toward the floor, raise the ribbed collar. Once the ribbed collar is at the threads, rotate it counterclockwise until it clicks over the lip on the Instrument Body and no gap remains between the ribbed collar and Instrument Body.
  - e. If placing the instrument in the Instrument Case, do not attach the Universal Implant Inserter Handle, as it has a separate set of brackets than the Implant Inserter.

### Sterilization:

The FORZA PTC Spacer System implants are supplied STERILE. They are sterilized using gamma irradiation sterilization. Do not re-sterilize.

### Sterilization in Orthofix Cases with Blue Wrap:

FORZA PEEK Spacer System instruments are compatible for use with the FORZA PTC Spacer System implants and are supplied NON-STERILE. Prior to use, all instruments should be placed in the appropriate Orthofix case which will be wrapped in a FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Preconditioning: 4 pulses
Drying time: 30 minutes		Exposure time: 4 minutes
Double wrapped		Drying time: 30 minutes
		Double wrapped

### Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer's instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix case (recommended 23¼" long x 11¼" wide container). Based on the rigid sterilization container size, the FORZA Discectomy Case may require removal of the lid prior to insertion into the sterilization container. The following sterilization cycle was validated:

Method: Steam
Cycle: Prevac
Temperature: 270°F (132°C)
Preconditioning: 4 pulses
Exposure time: 4 minutes
Drying time: 30 minutes

**Note:** Rigid sterilization containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10<sup>-6</sup>.

### Packaging:

The FORZA PTC Spacer System implants are sterilized using Gamma irradiation sterilization. Do not re-sterilize.

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

FORZA Spacer System instruments are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

### Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).

### Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

### Latex Information:

The implants, instruments and/or packaging material for the FORZA PTC Spacer System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

<b>Rx Only</b> Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
 See Instructions for Use	 Use By Date
 Orthofix.com/IFU	 Manufacturer
 Single Use Only Do Not Reuse	 Authorized Representative
 Catalogue Number	 Do Not Resterilize
 Provided Non-Sterile	 Serial Number
 Sterilized Using Irradiation	 Lot Number

Orthofix Inc.  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

Australian Sponsor  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

 Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

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