# INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

**Rx Only**  
CE 0086

 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

## Device System Name: **SKYHAWK<sup>®</sup> Lateral Plate System**

 Click directory below for desired language

<b>English</b> 	<b>2-3</b>
<b>Español</b> 	<b>4-5</b>

**Orthofix Inc.**  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

**Australian Sponsor**  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

**EC REP** Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

## English EN

### Device System Name: **SKYHAWK® Lateral Plate System**

**Description:**  
 The SKYHAWK Lateral Plate System consists of an assortment of non-sterile, single use, titanium alloy (Ti6Al4V ELI per ASTM F136) rigid plates and bone screws. The plates attach by means of screws to the anterolateral or lateral portion of the vertebral body of the thoracolumbar spine (T1-L5). The system includes instrumentation which assists in the surgical implantation of the devices. The system implants and instruments are provided non-sterile. They require sterilization prior to use.

**Indications for Use:**  
 The SKYHAWK Lateral Plate System is intended to be used as a non-pedicle lateral or anterolateral fixation system in skeletally mature patients and is intended to provide immobilization and stabilization of spinal segments as an adjunct to fusion in the treatment of the following acute and chronic instabilities of the thoracic and lumbar spine. It may be used from levels T1 to L5 with the following indications:

1. Degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies).
2. Spondylolisthesis.
3. Spinal stenosis.
4. Spinal deformities (i.e., scoliosis, kyphosis, and/or lordosis).
5. Tumor.
6. Pseudoarthrosis.
7. Failed previous fusion.
8. Trauma (i.e., fracture or dislocation).

**Contraindications:**  
 The SKYHAWK Lateral Plate System is contraindicated in patients with a systemic infection, with a local inflammation at the bone site, or with rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis. Do not use this system in patients with known or suspected metal allergies. Use of the system is also contraindicated in patients with any other medical, surgical or psychological condition that would preclude potential benefits of internal fixation surgery such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other disease, elevation of white blood cells or a marked shift in white blood cell differential count.

**Potential Adverse Events:**  
 Potential adverse events include, but are not limited to:

1. Early or late loosening of any or all of the components.
2. Disassembly, bending, and/or breakage of any or all of the components.
3. Foreign body (allergic) reaction to implants, debris, corrosion products, graft material, including metallosis, straining, tumor formation, and/or auto-immune disease.
4. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
5. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
6. Infection.
7. Vertebral body fracture at, above, or below the level of surgery.
8. Loss of neurological function, including paralysis (complete or incomplete).
9. Non-union, delayed union.
10. Pain, discomfort, or abnormal sensations due to the presence of the device.
11. Hemorrhage.
12. Cessation of any potential growth of the operated portion of the spine.
13. Death.

**Note:** Additional surgery may be necessary to correct some of these anticipated adverse events.

**Warnings and Precautions:**  
 The surgeon should be aware of the following when using implants:

1. The SKYHAWK Lateral Plate System is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.
2. Always orient the plate along the midline of the spine.
3. To optimize bony union, perform an anterior microdiscectomy or corpectomy as indicated.

4. To facilitate fusion, a sufficient quantity of autologous bone should be used.
5. Excessive torque applied to the screws when seating the plate may strip the threads in the bone.
6. Failure to achieve arthrodesis will result in eventual loosening and failure of the device construct.
7. **DO NOT USE DAMAGED IMPLANTS.**
8. Non-sterile; the system implants and instruments are provided non-sterile, and therefore must be sterilized before each use.
9. The health care provider must thoroughly clean the instruments prior to steam sterilization.
10. Single Use Only – system implants are SINGLE USE ONLY. No surgical implants should be reused. Any implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns that may lead to fatigue failure. Reuse could result in injury or require reoperation due to breakage or infection.
11. Do not combine system implantable components with those from any other system or manufacturer.
12. Do not combine implants of dissimilar metals, as it may result in galvanic corrosion.

**MRI Compatibility Information:**  
 The SKYHAWK Lateral Plate System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

**Cleaning:**  
 Implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluid it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned prior to use and after each re-use of the instruments. Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

None of the instruments in the system require disassembly prior to cleaning.

**From Point of Use:**  
 Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

**Preparation for Cleaning:**

1. Any instruments with moving parts (e.g., knobs, triggers and hinges) should be placed in the open position to allow better access of the cleaning fluid to the difficult to clean areas.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process. Use a soft cloth or plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning.

**Manual Cleaning:**

1. Because none of the instruments can be disassembled, use a soft cloth and/or a soft plastic bristle brush to remove any visible soil from the outside and inside of the instrument.
2. Prepare Vesphene® Ilse agent at the use-dilution recommended by label direction – 1 ounce per gallon or 1 mL of Vesphene® Ilse to each 128 mL of potable tap water per manufacturer's recommendations.
3. Bathe instruments in prepared room temperature solution as recommended by the detergent manufacturer.
4. Manually agitate instruments in Vesphene® Ilse solution for 15 minutes.
5. Scrub instruments with a soft plastic bristle brush if visible soil is noted and use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil.
6. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
7. Hang dry the device.
8. Visually inspect the instruments for visible soil.
9. If visible soil is noted, repeat the steps listed above.

**Automated Cleaning:**

1. Use a soft cloth or soft plastic bristle brush to remove any visible soil from the outside and inside of the instruments. Use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil from the inner lumens prior to automated cleaning.
2. Set up the washer/disinfectant detergent dose as ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer recommendations.
3. Remove the instruments from the instrument case and place them into the washer baskets.
4. Orient instruments in the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
  - a. Pre-Wash 1: cold potable water, 2 minutes
  - b. Enzyme/Detergent treatment:
    1. Spray, 20 seconds
    2. Soak, 1 minute
    3. Rinse cold potable water, 15 seconds
    4. Rinse cold potable water, 15 seconds
  - c. Wash ≥ 65°C, 2 minutes using Endozime AW Plus®
  - d. Rinse 1: hot potable water, 15 seconds
  - e. Rinse 2: hot potable water, 15 seconds
  - f. Rinse 3: hot potable water, 15 seconds
  - g. Rinse 4: hot potable water, 15 seconds
  - h. Thermal rinse ≥ 93°C, 1 minute
  - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
  - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
  - k. Dryat 115°C, 7 minutes
6. Visually inspect the instruments for visible soil
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

**Note:** Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

**Note:** Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix if instruments are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

**Sterilization:**

The SKYHAWK Lateral Plate System implants and instruments are supplied NON-STERILE. Prior to use, all implants and instruments should be placed in the instrument/implant case which will be either wrapped in an FDA cleared sterilization wrap or placed in a rigid sterilization container and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

**Sterilization in Orthofix Cases with Blue Wrap:**

- |                                |     |                                |
|--------------------------------|-----|--------------------------------|
| Method: Steam                  | or: | Method: Steam                  |
| Cycle: Gravity                 |     | Cycle: Prevac                  |
| Temperature: 270°F (132°C)     |     | Temperature: 270°F (132°C)     |
| Exposure time: 15 minutes      |     | Exposure time: 4 minutes       |
| Drying time: 30 minutes        |     | Drying time: 30 minutes        |
| Double wrapped                 |     | Double wrapped                 |
| (FDA cleared wrap recommended) |     | (FDA cleared wrap recommended) |

**Sterilization in Rigid Sterilization Containers:**

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix case (recommended 23¼" long x 11¼" wide container). The following sterilization cycle has been validated:

- Method: Steam
- Cycle: Prevac
- Temperature: 270°F (132°C)
- Exposure time: 4 minutes
- Drying time: 30 minutes

**Note:** Rigid sterilization containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10<sup>-6</sup>.

**Packaging:**

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The SKYHAWK Lateral Plate System instruments and implants are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

**Product Complaints:**

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness, and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).


**Further information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the SKYHAWK Lateral Plate System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

<b>Rx Only</b>	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use	 Manufacturer
	Orthofix.com/IFU	 Authorized Representative
	Single Use Only Do Not Reuse	 Serial Number
	Catalogue Number	 Lot Number
	Provided Non-Sterile	

 Orthofix Inc.  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

Australian Sponsor  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

 Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

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