 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**TDX™ Posterior Dynamic Stabilization System**

 Click directory below for desired language

<b>English</b> 	<b>2-3</b>
<b>Español</b> 	<b>4-5</b>
<b>Italiano</b> 	<b>6-7</b>
<b>Deutsch</b> 	<b>8-9</b>

**Orthofix Inc.**  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

**Australian Sponsor**  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

**EC REP** Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

## English EN

### Device System Name:

# TDX™ Posterior Dynamic Stabilization System

### Description:

The TDX Posterior Dynamic Stabilization System is a STERILE, single-use only motion preserving dynamic rod component used in conjunction with an Orthofix posterior pedicle screw system that allows a surgeon to build a spinal implant construct. The dynamic rod attaches to the vertebral body by means of screws and hooks to the non-cervical spine and consists of a polyurethane core secured in an implant-grade titanium housing.

### Indications for Use:

When used with posterior pedicle screw fixation systems in skeletally mature patients, i.e., the Firebird Spinal Fixation System or Spinal Fixation System (SFS), the TDX Posterior Dynamic Stabilization System is intended to provide immobilization and stabilization of non-cervical spinal segments as an adjunct to fusion in the treatment of the following acute and chronic instabilities using autogenous graft only:

1. Degenerative spondylolisthesis with objective evidence of neurologic impairment.
2. Kyphosis.
3. Failed previous fusion (pseudoarthrosis).

When used as part of a pedicle screw implant system, the TDX Posterior Dynamic Stabilization System is intended only for grade 3 or 4 spondylolisthesis at the fifth lumbar – first sacral (L5-S1) vertebral joint.

Please refer to the Firebird Spinal Fixation System or Spinal Fixation System (SFS) Instructions for Use (IFU) for a complete list of Indications for Use.

**Note:** For all indications, autogenous bone graft must be used.

### Contraindications:

Contraindications include, but are not limited to:

1. Morbid obesity.
2. Mental Illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia.
7. Patients unwilling or unable to follow post-operative care instructions.
8. If implanting the TDX Posterior Dynamic Stabilization System, any known allergy to titanium, polyurethane, or ethylene oxide residuals.
9. Any circumstances not listed under the heading indications.

### Potential Adverse Events:

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Device component fracture.
2. Loss of fixation.
3. Non-union.
4. Fracture of the vertebra.
5. Neurological injury.
6. Vascular or visceral injury.
7. Early or late loosening of any or all of the components.
8. Disassembly and/or bending of any or all components.
9. Foreign body (allergic) reaction to implants, debris, corrosion products, graft material, including metallosis, straining, tumor formation, and/or auto-immune disease.
10. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
11. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
12. Infection.
13. Pain, discomfort, or abnormal sensations due to the presence of the device.
14. Hemorrhage.
15. Cessation of any potential growth of the operated portion of the spine.
16. Death.

**Note:** Potential risks identified with the use of the device system may require additional surgery.

### Warnings and Precautions:

1. The safety and effectiveness of pedicle screw systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are: significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other condition are unknown.
2. The benefit of spinal fusions utilizing any pedicle screw fixation system has not been adequately established in patients with stable spines.
3. Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, loss of fixation, nonunion, fracture of the vertebra, neurological injury, and vascular or visceral injury.
4. Single use only.
5. The system instruments are sold non-sterile, and therefore must be sterilized before use.
6. The TDX dynamic rod is sold STERILE, and therefore should not be re-sterilized.
7. The system requires the use of autogenous bone graft.
8. Failure to achieve arthrodesis will result in eventual loosening and failure of the device construct.
9. DO NOT REUSE IMPLANTS. Discard used, damaged, or otherwise suspect implants.
10. The implantation of the system should be performed only by experienced spinal surgeons with specific training in the use of the dynamic rod as part of a pedicle screw spinal system.
11. Physicians/Surgeons should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact the performance of the system.
12. The safety and effectiveness of this device for the indication of spinal stabilization without fusion have not been established.
13. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.
14. The correct handling of the implant is extremely important. Implants should not be excessively or repeatedly bent, notched or scratched. These operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.

### MRI Compatibility Information:

The TDX Posterior Dynamic Stabilization System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

### Cleaning:

All instruments must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

The TDX Posterior Dynamic Stabilization System implants are provided STERILE and do not require cleaning.

### Sterilization:

The TDX Posterior Dynamic Stabilization System Instruments should be sterilized by the hospital using one of the following recommended cycles:

Method: Steam Cycle: Gravity Temperature: 250°F (121°C) Exposure time: 30 minutes	or:	Method: Steam Cycle: Prevac Temperature: 270°F (132°C) Exposure time: 8 minutes
--	-----	--

The TDX Posterior Dynamic Stabilization System implants are provided STERILE. They are sterilized using Ethylene Oxide sterilization. Do not re-sterilize.

### Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The TDX Posterior Dynamic Stabilization System instruments are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

TDX Posterior Dynamic Stabilization System implants are provided STERILE. Do not use if the package is opened or damaged, or if the expiration date has passed.

**Product Complaints:**

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).




**Further Information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the TDX Posterior Dynamic Stabilization System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** The TDX Posterior Dynamic Stabilization System has not been cleared by the US Food and Drug Administration and is not available in the United States.

	See Instructions for Use		Use By Date
	<a href="http://Orthofix.com/IFU">Orthofix.com/IFU</a>		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Catalogue Number		Do Not Resterilize
	Sterilized Using Ethylene Oxide		Serial Number
	Lot Number		

Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

**EC REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

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