 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**Ascent<sup>®</sup> Posterior Occipital  
Cervico-Thoracic (POCT) System**  
**Ascent<sup>®</sup> LE Posterior Occipital  
Cervico-Thoracic (POCT) System**

Click directory below for desired language

<b>English</b> 	<b>2-3</b>
<b>Español</b> 	<b>4-5</b>
<b>Français</b> 	<b>6-7</b>
<b>Italiano</b> 	<b>8-9</b>
<b>Deutsch</b> 	<b>10-11</b>
<b>Português</b> 	<b>12-13</b>

# INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

**Rx Only**  
CE 0086

Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

**EC/REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

## English EN

### Device System Name:

# Ascent® Posterior Occipital Cervico-Thoracic (POCT) System

# Ascent® LE Posterior Occipital Cervico-Thoracic (POCT) System

### Description:

The Ascent POCT and Ascent LE POCT Systems are temporary, multiple component systems comprised of a variety of non-sterile, single use components made of titanium alloy or cobalt chrome alloy that allow the surgeon to build a spinal implant construct. The Ascent POCT and Ascent LE POCT Systems consist of an assortment of rods, set screws, cross connectors, axial connectors, lateral offset adapters, multi-axial screws, hooks, plates, bone screws, and Songer Cables.

### Levels of Use:

When used in the occipito-cervico-thoracic spine, the Ascent POCT and Ascent LE POCT Systems may be used from the occiput to T3.

### Indications for Use:

When intended to promote fusion of the cervical spine and occipito-cervico-thoracic junction (occiput-T3), the Ascent POCT and Ascent LE POCT Systems are indicated for:

1. Degenerative disc disease (defined as neck pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies).
2. Spondylolisthesis.
3. Fracture/dislocation.
4. Spinal stenosis.
5. Atlanto-axial fracture with instability.
6. Occipito-cervical dislocation.
7. Tumors.
8. Revision of previous cervical spine surgery.

The occipital bone screws are limited to occipital fixation only. The use of the multi-axial screws is limited to placement in the upper thoracic spine (T1-T3) for the treatment of thoracic conditions only. They are not intended to be placed in the cervical spine. The lateral offset adapter is indicated for use in the upper thoracic spine (T1-T3). The hooks are intended to be placed from C1 to T3. The Songer Cable (titanium) System to be used with the Ascent POCT and Ascent LE POCT Systems allows for wire/cable attachment to the posterior cervical spine.

The Ascent POCT and Ascent LE POCT Systems can also be linked to the Orthofix Spinal Fixation System using the Axial or Parallel Rod Connector.

### Contraindications:

1. Morbid obesity.
2. Mental illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia.
7. Patients unwilling or unable to follow post-operative care instructions.
8. Any circumstances not listed under the heading Indications.

### Potential Adverse Events:

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Device component fracture.
2. Loss of fixation.
3. Non-union.
4. Fracture of the vertebra.
5. Neurological injury.
6. Vascular or visceral injury.
7. Early or late loosening of any or all of the components.
8. Disassembly and/or bending of any or all of the components.

9. Foreign body (allergic) reaction to implants, debris, corrosion products, graft material, including metallosis, straining, tumor formation, and/or autoimmune disease.
10. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation and/or pain.
11. Post-operative change in spinal curvature, loss of correction, height and/or reduction.
12. Infection.
13. Pain, discomfort or abnormal sensations due to the presence of the device.
14. Hemorrhage.
15. Cessation of any potential growth of the operated portion of the spine.
16. Death.

Note: Potential risks identified with the use of the device system may require additional surgery.

### Warnings and Precautions:

1. Single use only.
2. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.
3. Non-sterile; the screws, plates, rods, connectors, adapters, and instruments are sold non-sterile and therefore must be sterilized before each use.
4. Do not reuse implants; discard used, damaged or otherwise suspect implants.
5. Mixing of dissimilar metals can accelerate the corrosion process. Do not use the components of this system with implants of other material composition or components from different manufacturers unless specifically stated.

### MRI Compatibility Information:

The Ascent POCT and Ascent LE POCT Systems have not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The systems have not been tested for heating or migration in the MR environment.

### Cleaning:

All instruments and implants must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments and implants that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

### Sterilization:

#### Sterilization in Ascent Cases:

The Ascent POCT System should be sterilized by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

#### Sterilization in Ascent LE Cases:

The Ascent and Ascent LE POCT Systems' implants and instruments are supplied NON-STERILE. Prior to use, all implants and instruments should be placed in the instrumentation/implant case which will be wrapped in an FDA approved sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam (wrapped)	or:	Method: Steam (wrapped)
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Exposure time: 4 minutes
Drying time: 30 minutes		Drying time: 30 minutes

Ascent LE cobalt chrome occipital rods can be sterilized by placing them in the bottom of the Auxiliary Case. Ascent LE cobalt chrome straight rods can be sterilized by placing them in the bottom of the middle section of the Implant Case.

### Physician Information:

#### Patient Selection:

Patient selection is an extremely important factor in the success of implant procedures. It is important that the candidates be carefully screened and the optimal therapy selected.

**Preoperative:**

1. Carefully screen the patient, choosing only those that fit the indications described above.
2. Care should be exercised in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Store away from corrosive environments.
3. The construct should be assembled prior to surgery. An adequate inventory of implants of those expected to be used should be available at surgery.
4. All components and instruments should be cleaned and sterilized prior to use. Additional sterile components should be available in case of an unexpected need.

**Intraoperative:**

1. Instructions should be carefully followed.
2. Extreme caution should be used around the spinal cord and nerve roots.
3. The implant surface should not be scratched or notched since such actions may reduce the functional strength of the construct.
4. Bone cement should not be used, as it will make removal of the components difficult or impossible.
5. Before closing soft tissue, check each screw to make sure that none have loosened.

**Postoperative:**

1. Detailed instructions should be given to the patient regarding care and limitations, if any.
2. To achieve maximum results, the patient should not be exposed to excessive mechanical vibration. The patient should not smoke or consume alcohol during the healing process.
3. The patient should be advised of their limitations and taught to compensate for this permanent physical restriction in body motion.
4. If a non-union develops, or if the components loosen, the devices should be revised or removed before serious injury occurs. Failure to immobilize the non-union, or a delay in such will result in excessive and repeated stresses on the implant. It is important that immobilization of the spinal segment be maintained until fusion has occurred.
5. The implants are temporary internal fixation devices. Internal fixation devices are designed to stabilize the spine during the normal healing process. After the spine is fused, the devices serve no functional purpose and should be removed.

**Patient Information:**

The temporary internal fixation devices used in your recent spinal injury are metallic implants that attach to the bone and aid in the healing of bone grafts. These implants have been shown to be valuable aids to surgeons in the treatment of bony fusions. These devices do not have the capabilities of living bone. Intact living bone is self-repairing, flexible and occasionally breaks and/or degrades. The anatomy of the human body places a size limitation on any artificial fixation device used in surgery. The maximum size limitation increases the chances of the mechanical complications of loosening, bending, or breaking of the devices. Any of these complications could result in the need for additional surgery. Accordingly, it is very important that you follow the recommendations of your physician. Use braces as instructed. By following these instructions, you can increase your chances of a successful result and reduce your risk of injury and/or additional surgery.

**Packaging:**

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The Ascent POCT and Ascent LE POCT Systems' instruments and implants are provided in modular cases specifically intended to contain and organize each system's components. The systems' instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

**Product Complaints:**

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).


**Further Information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the Ascent and Ascent LE POCT Systems are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician

<b>Rx Only</b>		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Catalogue Number
	Orthofix.com/IFU		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Provided Non-Sterile		Serial Number
	Lot Number		

 Orthofix Inc.  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

Australian Sponsor  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

 Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

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