



Innovation takes root in great standards

Chimaera Hip Fracture System™
designed to leverage decades of knowledge
and experience in hip nailing



Indicated for treatment of **stable and unstable** pertrochanteric, intertrochanteric and subtrochanteric **fractures**


All-in-one instrument tray
Intuitive **color coded** instrument selection

Instrumentation allows reproducible, easy and **accurate measurement** and may **reduce related x-ray exposure**

Supplementary lag screw **anti-rotational stability**

Self-locking sliding lag screw

Designed to **reduce lateral pain** during weight-bearing and the related **reoperation rate**



1:1 scale

TECHNICAL DETAILS	
Short nails length:	180mm
Long nails length:	280-460mm (20mm increments)
Angulations:	125°, 130°
Metaphyseal angle:	5°
Proximal diameter:	15.5mm
Distal diameter:	10mm, 11mm
Cannulated guided end caps	

Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



www.orthofix.com

HF-1705-PL-E0 A 01/18