


PHOENIX[®] CDX[™]
Minimally Invasive Spinal Fixation System - Compression Distraction Extension



Advancing Phoenix[®] One Level at a Time

As an extension to the already versatile Phoenix Minimally Invasive System, Phoenix Compression Distraction Extension offers even greater implant and instrument options for complex MIS deformity spinal fixation.

Mono-Axial Phoenix Body

- Unique MIS modular screw design that can be used to facilitate direct vertebral rotation.

Mono-Axial Body





Fast Drive

- Rotates set screw driver approximately three times per turn



Rod Inserter

- New design that allows for increased visualization and improved functionality




DVR Tower

- Perform Direct Vertebral Rotation percutaneously with use of ratcheting connectors and rotation handles
- Perform final tightening with use of counter torque wrench
- Holes on shaft for compressor/distractor securement



MIS Counter Torque Tower

- Counter torque through opening of Phoenix Body without increasing incision diameter.
- Perform final tightening with use of counter torque wrench
- Holes on shaft for compressor/distractor securement



Compressor/Distractor

- Multi-Level parallel Compression and Distraction
- Longer Throw
- Interchangeable Tips

Length of Throw



60mm

60mm Tip

Compression: 15cm

Distraction: 13cm



30mm

30mm Tip

Compression: 12cm

Distraction: 10cm



Straight

Straight Tip

Compression: 9cm

Distraction: 7cm

Trays

Part #	Description	Qty
Top Level		
20-0320	Compression/Distractor Fixture	1
20-0380	Compressor	1
20-0381	Distractor	1
20-0382	Compression/Distractor Tip	2
20-0383	30mm Offset Compression/Distractor Tip	2
20-0384	60mm Offset Compression/Distractor Tip	2
20-7450	450mm Rod, CoCr	2


Implants, Order By Request

20-8090	Mono-Axial Open Body, Standard
20-4200	200mm Prelordosed Rod


Instruments, Order By Request

20-0300	Mono-Axial Screw Driver
20-0310	Fast Drive
20-0311	Handle Adapter
20-0313	Rod Inserter
20-0321	MIS Counter Torque Tower
20-0330	DVR Tower, Short - Standard
20-0340	Set Screw Driver
51-1986	Extra Large Ratcheting Connector
51-1486	Rotation Handle

Distributed by:

 Manufactured by:
Orthofix
3451 Plano Parkway
Lewisville, Texas 75056-9453 USA
214-937-2000

 0086 **RX Only**



Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for full information on indications for use, contraindications, warnings, precautions, adverse reactions information and sterilization.

1.888.298.5700
www.orthofix.com


ORTHOFIX®