




**SKYHAWK™**

Lateral Interbody Fusion System  
Lateral Plate System





## Complete Solution for Lateral Fusion

Lateral Lumbar Interbody Fusion is an approach to spinal fusion in which the surgeon accesses the intervertebral disc space using a surgical approach from the patient's side (lateral) rather than from the front (anterior) or the back (posterior). The lateral approach provides an alternate route to the spine that disturbs fewer structures and tissues. The SKYHAWK™ Lateral Interbody Fusion System and SKYHAWK Lateral Plate System provide a complete solution for the surgeon to perform a Lateral Lumbar Interbody Fusion (L2-L5).

### Lateral Interbody Fusion System

#### Implant

- Standard Interbody offerings of 18mm and 22mm widths (additional sizes available upon request)
- 0° and 8° standard implants (12° available upon request)
- Tantalum markers located at both ends and the middle of the implant for clear fluoroscopic visualization
- Large graft window for packing bone grafting material



### Lateral Plate System

#### Implant

- Low profile Lateral Plate
- Set screws seat flush in the plate
- Plate sizes coordinate with Interbody heights



#### Screw

- 5.5mm Primary Screw: 30mm-60mm (5mm increments)
- 6.5mm Rescue Screw: 30mm-60mm (5mm increments)





# QUADR<sup>x</sup>

LATERAL RETRACTOR SYSTEM

## Simple Mechanical Functions Make QUADR<sup>x</sup>™ a Seamless Addition to any Lateral Procedure

- Implements a straightforward 4 blade design
- Ergonomic rectangular aperture opening provides adaptability to both patient anatomy and surgeon technique
- Innovative blade design requires no sequential dilation prior to retractor placement over the guide wire





## Instruments

Part #	Description	Qty Set
<b>Retractor Case</b>		
71-0110	Retractor Case	1
71-1000	Retractor	1
71-1100	Blade, 100mm	5
71-1110	Blade, 110mm	5
71-1120	Blade, 120mm	5
71-1130	Blade, 130mm	5
71-1140	Blade, 140mm	5
71-1150	Blade, 150mm	5
71-1160	Blade, 160mm	5
71-1500	Retractor Driver	1
71-4100	Intradiscal Shim	2
71-4200	Bone Pin	2
71-4300	Shim Insertor	1

### Accessory Case

71-0120	Retractor Accessory Case	1
71-2000	Table Arm	1
71-2010	Table Clamp	1
71-2030	Retractor Light Bracket	1
71-3030	Fiberoptic Light Cable	1
71-5010	Fluoroscopic Targeting Crosshair	1
71-5015	Guide-Wire	10
71-5040	Dilator Holder	1

### Sterile Packaged


71-3040	Fiberoptic Light Source	1
7332-99	Neuromonitoring Pack	1

## Implants

Part #	Description	Qty Set
<b>Lateral Plate Case</b>		
24-2001	Set Screw	8
24-2006	6mm Plate	4
24-2008	8mm Plate	4
24-2010	10mm Plate	4
24-2012	12mm Plate	2
24-2014	14mm Plate	2
24-5530	5.5mm x 30mm Bone Screw	4
24-5535	5.5mm x 35mm Bone Screw	6
24-5540	5.5mm x 40mm Bone Screw	6
24-5545	5.5mm x 45mm Bone Screw	6
24-5550	5.5mm x 50mm Bone Screw	6
24-5555	5.5mm x 55mm Bone Screw	6
24-5560	5.5mm x 60mm Bone Screw	4
24-5630	6.5mm x 30mm Bone Screw	4
24-5635	6.5mm x 35mm Bone Screw	6
24-5640	6.5mm x 40mm Bone Screw	6
24-5645	6.5mm x 45mm Bone Screw	6
24-5650	6.5mm x 50mm Bone Screw	6
24-5655	6.5mm x 55mm Bone Screw	6
24-5660	6.5mm x 60mm Bone Screw	4

### Interbody 18W Implant Case

27-2208	18W x 40L x 0 Deg x 8H Interbody	2
27-2210	18W x 40L x 0 DEG x 10H Interbody	2
27-2212	18W x 40L x 0 DEG x 12H Interbody	2
27-2214	18W x 40L x 0 DEG x 14H Interbody	2
27-2238	18W x 40L x 8 DEG x 8H Interbody	2




Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Proper surgical procedure is the responsibility of the medical professional. Please refer to the "Instructions for Use" supplied with the product for full information on indications for use, contraindications, warnings, precautions, adverse reactions information and sterilization.

1.888.298.5700  
[www.orthofix.com](http://www.orthofix.com)

SH-1405-PL-US © Orthofix Holdings, Inc. 12/2014

 Orthofix  
 3451 Plano Parkway  
 Lewisville, Texas 75056 U.S.A.

## Implants

Part #	Description	Qty Set
<b>Interbody 18W Implant Case (cont)</b>		
27-2240	18W x 40L x 8 DEG x 10H Interbody	2
27-2242	18W x 40L x 8 DEG x 12H Interbody	2
27-2244	18W x 40L x 8 DEG x 14H Interbody	2
27-2308	18W x 45L x 0 DEG x 8H Interbody	2
27-2310	18W x 45L x 0 DEG x 10H Interbody	2
27-2312	18W x 45L x 0 DEG x 12H Interbody	2
27-2314	18W x 45L x 0 DEG x 14H Interbody	2
27-2338	18W x 45L x 8 DEG x 8H Interbody	2
27-2340	18W x 45L x 8 DEG x 10H Interbody	2
27-2342	18W x 45L x 8 DEG x 12H Interbody	2
27-2344	18W x 45L x 8 DEG x 14H Interbody	2
27-2408	18W x 50L x 0 DEG x 8H Interbody	2
27-2410	18W x 50L x 0 DEG x 10H Interbody	2
27-2412	18W x 50L x 0 DEG x 12H Interbody	2
27-2414	18W x 50L x 0 DEG x 14H Interbody	2
27-2438	18W x 50L x 8 DEG x 8H Interbody	2
27-2440	18W x 50L x 8 DEG x 10H Interbody	2
27-2442	18W x 50L x 8 DEG x 12H Interbody	2
27-2444	18W x 50L x 8 DEG x 14H Interbody	2
27-2508	18W x 55L x 0 DEG x 8H Interbody	2
27-2510	18W x 55L x 0 DEG x 10H Interbody	2
27-2512	18W x 55L x 0 DEG x 12H Interbody	2
27-2514	18W x 55L x 0 DEG x 14H Interbody	2
27-2538	18W x 55L x 8 DEG x 8H Interbody	2
27-2540	18W x 55L x 8 DEG x 10H Interbody	2
27-2542	18W x 55L x 8 DEG x 12H Interbody	2
27-2544	18W x 55L x 8 DEG x 14H Interbody	2

### Interbody 22W Implant Case

27-3208	22W x 40L x 0 DEG x 8H Interbody	2
27-3210	22W x 40L x 0 DEG x 10H Interbody	2
27-3212	22W x 40L x 0 DEG x 12H Interbody	2
27-3214	22W x 40L x 0 DEG x 14H Interbody	2
27-3238	22W x 40L x 8 DEG x 8H Interbody	2
27-3240	22W x 40L x 8 DEG x 10H Interbody	2
27-3242	22W x 40L x 8 DEG x 12H Interbody	2
27-3244	22W x 40L x 8 DEG x 14H Interbody	2
27-3308	22W x 45L x 0 DEG x 8H Interbody	2
27-3310	22W x 45L x 0 DEG x 10H Interbody	2
27-3312	22W x 45L x 0 DEG x 12H Interbody	2
27-3314	22W x 45L x 0 DEG x 14H Interbody	2
27-3338	22W x 45L x 8 DEG x 8H Interbody	2
27-3340	22W x 45L x 8 DEG x 10H Interbody	2
27-3342	22W x 45L x 8 DEG x 12H Interbody	2
27-3344	22W x 45L x 8 DEG x 14H Interbody	2
27-3408	22W x 50L x 0 DEG x 8H Interbody	2
27-3410	22W x 50L x 0 DEG x 10H Interbody	2
27-3412	22W x 50L x 0 DEG x 12H Interbody	2
27-3414	22W x 50L x 0 DEG x 14H Interbody	2
27-3438	22W x 50L x 8 DEG x 8H Interbody	2
27-3440	22W x 50L x 8 DEG x 10H Interbody	2
27-3442	22W x 50L x 8 DEG x 12H Interbody	2
27-3444	22W x 50L x 8 DEG x 14H Interbody	2
27-3508	22W x 55L x 0 DEG x 8H Interbody	2
27-3510	22W x 55L x 0 DEG x 10H Interbody	2
27-3512	22W x 55L x 0 DEG x 12H Interbody	2
27-3514	22W x 55L x 0 DEG x 14H Interbody	2
27-3538	22W x 55L x 8 DEG x 8H Interbody	2
27-3540	22W x 55L x 8 DEG x 10H Interbody	2
27-3542	22W x 55L x 8 DEG x 12H Interbody	2
27-3544	22W x 55L x 8 DEG x 14H Interbody	2

  
**ORTHOFIX®**